

Acutis Reveal™ UTI collection devices

To ensure your practice's Acutis Reveal™ UTI samples are processed in the most efficient way possible, please utilize the following guidelines when selecting the proper collection device for each stage of the UTI test.

For confirmative
Acutis Reveal™ UTI PCR testing

BD Vacutainer® Plus C&S preservative tube 4.0 mL




- Tube dimensions:
 - Plus plastic tube
 - Size: 13 x 75 mm
 - Conventional stopper
- 4.0 mL
- Boric Acid
Sodium Formate
Sodium Borate preservative

Ref# 364951

For urinalysis screening

BD Vacutainer® Plus urinalysis tube



- Tube dimensions:
 - Plus plastic tube
 - Size: 16 x 100 mm
 - Conventional stopper
- 8.0 mL
- Ethyl Paraben
Sodium Propionate
Chlohexidine preservative

Ref# 364992